




CQ-00971-D2O-DEPT135-271009-298 K-500 MHz  
 DEPT135-mmmod-clem\_D2O /opt/topspin ravi {1 C12 - 136}



NAME CQ-00971-D2O-271009  
 EXPNO 5  
 PROCNO 1  
 Date\_ 20091031  
 Time 15.43  
 INSTRUM spect  
 PROBHD 5 mm CPTXO 13C  
 PULPROG dept135  
 TD 65536  
 SOLVENT D2O  
 NS 256  
 DS 16  
 SWH 29761.904 Hz  
 FIDRES 0.454131 Hz  
 AQ 1.1010548 sec  
 RG 203  
 DW 16.800 usec  
 DE 30.00 usec  
 TE 297.9 K  
 CNST2 145.0000000  
 D1 2.00000000 sec  
 D2 0.00344828 sec  
 D12 0.00002000 sec  
 TD0 1

===== CHANNEL f1 =====  
 NUC1 13C  
 P1 10.00 usec  
 P2 20.00 usec  
 PL1 3.00 dB  
 PL1W 36.25705338 W  
 SFO1 125.6974365 MHz

===== CHANNEL f2 =====  
 CPDPRG2 waltz16  
 NUC2 1H  
 P3 15.70 usec  
 P4 31.40 usec  
 PCPD2 100.00 usec  
 PL2 1.25 dB  
 PL12 17.33 dB  
 PL2W 19.96668053 W  
 PL12W 0.49238619 W  
 SFO2 499.8419994 MHz  
 SI 32768  
 SF 125.6848680 MHz  
 WDW EM  
 SSB 0  
 LB 1.00 Hz  
 GB 0  
 PC 1.40